

4C Lab – Optical Instruments

Goal: To examine the properties of thin lenses and construct optical instruments

Equipment List:

- Optics bench
- Laser bench (optional)
- Box of lenses
- Component holders (3)
- Light box with target
- 3 x 5 card
- 2 meter stick

Pre-lab exercise: Determine the position of the image formed with the telescope

Background: The thin lens equation:


$$\frac{1}{s} + \frac{1}{s'} = \frac{1}{f}$$

Describes the relationship between the image distance, s , the object distance s' and the focal length f . This equation also holds for a system of lenses using the image formed by the first as the object for the second.

Procedure – Part 1: The Simple Magnifier

1. Find the 48mm lens. Measure the focal length by forming the real image of a distant, real object (a tree outside, the moon, your lab partner).
2. Using the crossed arrow target on the light box, find the focal length using the thin lens equation. Measure s , s' . Find f .
3. Measure y and y' . Find the lateral magnification, $m = -s'/s = y'/y$ using both ratios and compare.
4. Repeat for the 252 mm lens. This one is a bit trickier, but with a bit of thought and patience completely doable. (Hint: Make sure the lens is not closer than 252 mm to the target and start with the 3 x 5 card 1.5 m away)


Part 2: A compound microscope – Looking at small things at short distances



- just
- Place the 48 mm lens, the objective, on the track so that the crossed arrow target is ^{just} beyond the first focal point. Find the image. It will be real, enlarged and inverted.
 - Place the 127mm lens, the eyepiece, so that its first focal point coincides with the image formed by the 48mm lens.
 - Determine the tube length of your microscope L . This is the distance between the second focal point of the objective and the first focal point of the eyepiece.
 - Try to view the crossed arrow target through the eyepiece. Measure how far away the lens is from your eye.
 - Find the magnifying power of the microscope: The lateral magnification of the objective is $m_o = y'/y = -L/f_o$. The angular magnification of the eyepiece is $M_e = x_{np}/f_e$, where x_{np} is the near point of the viewer (the nearest point the viewer can focus), and f_e is the focal length of the eyepiece. The magnifying power of the microscope is the product of the the two:

$$M = m_o M_e = -\frac{L}{f_o} \frac{x_{np}}{f_e}$$

Part 3. The Telescope – Looking at large distant objects.



- Place the 252 mm lens, the objective, on the bench.
- Looking at a distant object with the lens (something outside, the open door, the exit sign?) find the position of image using the 3 x 5 card. Measure the height of the object on the 3 x 5 card. This is y' .
- Remove the 3 x 5 card and put the 48 mm lens at one end. The distance between the two lenses should be exactly the sum of the focal lengths.
- Try to view something using this telescope. It is challenging to look through both lenses at once, but give it a try. Choosing an object far away helps. Is the image upright or inverted?
- Calculate the angular magnification of the objective:

$$\tan \theta_o = -\frac{y'}{f_o} \approx \theta_o$$

and the angular magnification of the eyepiece:

$$\tan \theta_e = -\frac{y'}{f_e} \approx \theta_e$$

- The magnification of the telescope is:

$$M = \frac{\theta_e}{\theta_o} = -\frac{f_o}{f_e}$$

- Estimate the magnification you see and compare with the calculated values.